

torneo® II – active hoist

- Optimal kinematics – lifts and lowers patients gently
- Transfers safely into a standing as well as a seating position
- Patient's remaining mobility is optimally preserved
- Relieves back strain for caregivers
- Can be used optimally in small rooms and rooms with narrow doors
- Ergonomic active hoist with low total weight for user-friendly operation
- Highly stable, a foot lever mechanically inclines the carrier
- Maximum load: 150 kg

Active hoist
torneo® II



How the
torneo® II works:
www.youtube.de

foldy® – the original

- Easy, uncomplicated handling
- Easy to fold and transport, no tools necessary
- Lifts patients effortlessly, even from the floor
- The carrier is inclined via a foot lever (electronic option)
- Low ground clearance, suitable for aks-care beds
- Maximum load: 150 kg
- Also available in mini and micro versions (with a maximum load of 130 kg model), also as foldy® XL with a maximum load of 185 kg



clino® II – the quick

- Quick motor
- Easy running castors, easy to manoeuvre and push, even in small rooms
- Easy to fold and transport, no tools necessary
- Lifts patients effortlessly from the floor
- The carrier is inclined via a foot lever
- Maximum load: 150 kg
- Also available and stable for bigger and heavier patients: clino® XL up to 185 kg maximum load with strong and fast motor and goliath® up to 250 kg maximum load, ideal for care facilities

clino® II



Standard hoists

aks from Troisdorf, Germany

Our products are developed and produced in Troisdorf, Germany. We have been well-known manufacturers of technical care aids across Germany and Europe. Patients, caregivers and retailers value and approve of our products. Our patient hoists are absolutely safe, reliable and durable.

Medical aids for the care by aks:

- Mobile patient hoists
- Transfer aids
- Care beds
- Antidecubitus alternating pressure systems, foam mattresses and seat cushions
- Bariatric care products

aks

aks GmbH
D-53842 Troisdorf
T +49 (0) 2241/9474-0
www.aks.de

A company of the non-profit organisation EIN HERZ LACHT
www.einherzlacht.org

Advice and sale by:



Patient hoists

aks

We reserve the right to make changes and colour differences. Subject to printing errors. Stand 10/2018



aks active hoist

- torneo® II and torneo® compact
- dualo® mini, dualo®, dualo® XL with active hoist attachment

aks standard hoist

- foldy® micro, foldy® mini, foldy®, foldy® XL
- dualo® mini, dualo®, dualo® XL with standard hoist attachment
- clino® II, clino® XL
- goliath® for bariatric care

aks patient hoists

Back-friendly transfer and care systems

aks patient hoists

We offer patient hoists according to various requirements: patients of different sizes, different body weights and varying room sizes.

With mobile patient hoists, we differentiate between active and standard hoists.

Active hoists can be used to stand and transfer (in standing) patients with some remaining mobility. Slings offer support and stability during lifting and transfer.

Standard hoists are used, for example, to transfer patients from the bed to a chair or from the wheelchair to the toilet. The transfer is made using slings. The transfer is generally carried out while seated, but it is also possible to lift from the floor.

Our **patient sling program** provides a number of different slings for various pathologies, constitutions and care requirements. Slings for wet and dry areas, for simple transfers out of the bed or to use the toilet or bathtub.

dualo® – the versatile

- Versatile: Switch easily between the standard hoist and the active hoist as you wish, no tools necessary. Various optional and individual conversion kits.
- Multifunctional: efficient use of the standard and/or active hoist's various functions, e.g. lifting, straightening up and lowering patients.
- Economic and effective: low acquisition costs and broad range of utilisation in needs-orientated application as a dual lifter – particularly practical in care facilities.
- Maximum load: 150 kg



The dualo® is quickly modified



foldo® – the unique

- "Two-in-one concept" – **standard hoist and bathtub hoist all-in-one**
- Can be lowered to the bottom of the bathtub
- Lateral lifting for bathtub use
- No need for an additional bathtub, wall or ceiling hoist
- "Dock" the hoist conveniently without any structural modifications
- Bathtub support and telescopic jib provide bathroom safety by stable support at the bathtub edge.
- Maximum load: 130 kg



foldo® standard hoist with swivel arm and telescopic jib

Accessories for standard hoist



aks accessories for aks-standard hoist

- Hoist scales
- Removable battery pack and wall-loading stations
- Different spreader bars
- Wide patient sling-program for varying requirements in daily care

Transfer aids

The raising belt offers back-friendly support in uprighting and setting down patients. In doing so, the patient can hold on to the walking aid belt provided by the caregiver. Wearing the walking aid belt, the patient is stabilised by the caregiver when walking.

Slide sheet for lateral transfer reduces friction and can be used to move immobile patients without causing back strain.



Slide sheet for lateral transfer



The different types of hoists by aks:

- Active hoists
- Standard hoists

aks standard hoist with the PLUS of functions

- **dualo®** – from active to standard hoist, quickly modified with the optional conversation kit
- **foldo®** – from standard to a bathtub hoist with the swivel arm and the telescopic jib

Transfer aids

- Raising belt
- Walking belt
- Slide sheet for lateral transfer