

aks D4 low entry

The low entry bed for the prevention of falls



More and more people need safety and protection from falling out of bed due to the progression of their illness, such as dementia or Alzheimer's disease. The low entry height of only 22 cm contributes significantly to fall prevention, without having to subject the patient to measures that restrict his or her freedom. Likewise, the low entry bed offers smaller people independent mobility on a daily basis.

The advantages of the aks D4 low entry:

- Low positioning of persons, with an entry/exit height from as low as 22 cm
- The low height provides the patient with safety and a sense of well-being
- Suitable for patients who are afraid of falling out of the care bed (fall prevention)
- As an aid to mobility for small persons
- Maximum variability due to a four-section lying area made of metal mesh or wood sprung slats
- Maximum functionality: Height, head and foot parts electrically adjustable; the latter with knee bend is electrically adjustable from the bed by hand control without nursing assistance
- Twist-proof side railings integrated on both sides, easy assembly and dismantling by one person, particularly light head and foot sections
- Hand control with integrated lock function
- Transport and storage aid for unpackaged delivery
- Including uplifter and grab handle

- 3-year warranty on electrical components
- Castors with individual wheel locking
- Max. patient weight 150 kg, safe working load of 185 kg

Options and accessories

- Side bar height extender, e.g. for use with high cell mattresses over 14 cm or 16 cm, depending on the type of lying area
- Side bar padding
- Bed and mattress extension of 20 cm for the accommodation of larger patients
- Infusion holder
- Foot-end lowering (FEL)
- Unpackaged delivery on transport aid



Scan to see how the knee bend works.

 @aks_GmbH