







with a conversion kit, active and standard plugged together for storage



dualo® XL

strong – flexible – multifunctional

The well established, variably deployable dualo® patient lift has been designed for larger or heavier patients in an XL format. It is used for straightening and transferring a person in combination with aks hoist straps in a care setting at home or in nursing homes.

The advantages of the dualo® XL

- economy: lower acquisition costs
- effectiveness: high frequency of use, needs-oriented use as a combined hoist, especially in nursing homes
- the transfer/straightening of larger or heavier patients is also possible. A correspondingly powerful motor is required for safe lifting and lowering
- · max. load capacity 185 kg
- high stability as a result of the corresponding design
- multi-functionality: Optimum use of the versatile functions of a standard and/or straightening hoist, e.g. lifting, straightening and turning people. Support, also from the floor
- Flexibility: easy conversion of the active hoist (Fig.1) to a standard hoist (Fig.3) and vie versa. Simple plug-in function of the stainless steel rectangular tubes for guick, tool-free conversion. The conversion kit is optional
- the base of the chassis is hinged (Fig. 5)
- high-quality drive components from LINAK®
- new ergonomic handle
- high-quality, ball-bearing castors with threadguards and excellent rollability, easy manoeuvrability
- mechanical emergency lowering, emergency stop switch, acoustic warning against total discharge
- colour: light grey

Options and accessories

- additional battery pack, separate wall charger, calf strap, hoist scales, hoist straps
- optionally available with electrical spreader bars (dualo® e XL)

Order numbers

dualo® XL active (without strap) – 87311 dualo® e XL active (without strap) – 87331 dualo® XL standard (without frame and strap) – 87310 dualo® e XL standard (without frame and strap) – 87330

For further order numbers, see price list

