



aks sliding and turning sheet


The aks sliding and turning sheets are made of silicone-coated fabric that effectively reduces the friction resistance during the repositioning of patients. This makes the corresponding activities of care personnel easier. Different sizes of the aks sliding and turning sheets can be used depending on the restriction of the patient's mobility. The body regions to be supported are decisive here.

The large number of possible uses covers repositioning of normal, heavy and very heavy immobile patients from abdomen, back or side position to any required position. Good patient care simply becomes more feasible.

Technical Data

Outer material: silicone-coated fabric, green
washable at 60° C; do not dry in a dryer

Size 195 x 70 cm
Size 140 x 110 cm
Size 130 x 70 cm
Size 80 x 70 cm
Size 48 x 48 cm



aks sliding and turning sheet with grip loops

This design combines all positive product characteristics of the aks sliding and turning sheets and may potentially be used for particularly heavy and very heavy patients. The grip loops here provide as low friction resistance as possible while in use and thus make movement of these patients possible. With more and more people suffering from adiposity, the use of such sheets is becoming increasingly more important.

Technical Data

Outer material: silicone-coated fabric, green
Grip loops: Polyester sling, red
washable at 60 °C; do not dry in a dryer

Size 195 x 70 cm
Size 140 x 110 cm
Size 150 x 70 cm

