



## saniflow® II

### Large-cell, standard class anti-decubitus alternating pressure system

The standard class already starts with a large-cell mattress overlay system which provides the patient with pleasant relief using alternating pressure.

#### The advantages of the saniflow® II:

- Overlay mattress system
- Easy operation: continuously variable pressure setting according to the patient weight, stable folding metal spreader bar for attachment to the bed, simple filter replacement
- Visual display for reaching the required pressure
- Visual alarm in the event of pressure loss (e.g. leak)
- Highly elastic foam strips on the rear side and soft feet for noise insulation
- Cycle time 12 minutes
- Alternating two-chamber system for intermittent pressure relief
- 17/19 cells, can be replaced individually
- Cell height 13 cm
- Mattress cover made of nylon/PVC fastened to the mattress overlay at the head and foot ends with press studs, water-impermeable
- With fold-over ends for fastening to the foam mattress underlay

- Fast ventilation of the mattress using end plugs with pull ring
- Secure snap-in connection of the hoses on the pump unit
- Available in 90/100 cm x 200/220 cm

Patient weight 40-120 kg

Decubitus grade up to and incl. Grade II (according to EPUAP)