



saniflow® III

Large-cell, medium standard class anti-decubitus alternating pressure system with static function

This further developed version enables persons with up to 130 kg body weight and up to and including Decubitus Grade III according to EPUAP to use the advantages of a medium standard class alternating pressure system.

The advantages of the saniflow® III:

- Overlay mattress system
- With static function and thus suitable for patients in pain
- Easy operation: continuously variable pressure setting according to the patient weight; simple filter replacement
- Visual display for reaching the required pressure
- Visual and audible alarms in the event of pressure loss (e.g. leak)
- Switch for resetting the audible alarm
- Highly elastic foam strips on the rear side and soft feet for noise insulation
- Cycle time 12 minutes
- Alternating two-chamber system for intermittent pressure relief
- 17/19 cells, can be replaced individually

- Cell height 13 cm
- Mattress cover made of nylon/PVC fastened to the mattress with all-round zip; water-impermeable
- With diagonal rubber strips routed over the corners on the mattress underside for fixing to the foam mattress underlay
- Fast ventilation of the mattress using end plugs with pull ring
- Secure snap-in connection of the hoses on the pump unit
- Available in 90/100 cm x 200/220 cm

Patient weight 40-130 kg
Decubitus grade up to and incl. Grade III (according to EPUAP)