



saniflow® IV

Large-cell enhanced standard class anti-decubitus alternating pressure system with static and care function

As a mattress replacement system with a sophisticated range of functions, saniflow® IV meets the high requirements of the enhanced standard. Patients with decubitus up to and including Grade IV according to EPUAP receive the best standard care for decubitus therapy.

The advantages of the saniflow® IV:

- Mattress replacement system
- Care function enables a maximum filled static support surface with a time limit of 30 minutes and automatic changeover to the previously selected alternating pressure or static operating mode and the selected pressure level.
- 8-stage pressure setting according to the patient weight
- Visual display for reaching the required pressure
- Visual and audible alarms in the event of pressure loss (e.g. leak)
- Switch for resetting the audible alarm
- Cycle time can be set from 10-25 minutes
- Mattress replacement system with removable foam underlay (4 cm)
- Alternating two-chamber system for intermittent pressure relief
- 17/19 cells, can be replaced individually

- Cell height 13 cm
- Mattress cover made of nylon/PVC fastened to the mattress with all-round zip; water-impermeable
- With diagonal rubber strips routed over the corners on the mattress underside for fixing to the horizontal area
- Fast ventilation of the mattress using end plugs with pull ring
- Secure snap-in connection of the hoses on the pump unit
- Available in 90/100 cm x 200/220 cm

Patient weight 40-130 kg

Decubitus grade up to and incl. Grade IV (according to EPUAP)